



INNOVATIVE GLOVES

Our plant is a dream of professionals and perfectionist who have been involved in the latex dipping industry since the last 30 years. All of our products not only meet, but also exceed demanding international standards and regulations. To name a few we have surpassed USFDA requirements, EN 46002 and CE Mark Certification from TUV Germany.



Innovative Gloves is Thai owned company, manufacturing premium latex examination gloves located 2 miles from Hatyai International Airport in Thailand. Our facility has total area of 56,000 sq.m. with a covered area of over 10,000 sq.m.

We currently have a capacity to produce 30 million gloves a month and are expanding to 60 million gloves a month in the near future. We practically supply to customers from every corner of the world - in particular, the U.S, South America and Europe.



INNOVATIVE GLOVES

Salient Features of Our Plant and Products

- Our dipping lines are fully fabricated with stainless steel, to ensure that our products have negligible impurities, reduce risks of contamination and geared against corrosion.
- Our gloves are produced using Shinko formers, which are renowned in the industry for their comfort, superior and universal fittings.



- Our machines are designed 100 meters in length with 6 tiers to produce very high and consistent quality products with 3 pre-leaching and post leaching tanks of 6 meters each. These significantly reduce and extract maximum protein, which is the main cause of allergy.

- We have all 5 sizes i.e. XS, S, M, L and XL to cater for our customers' needs.
- We are capable of producing long cuff beaded gloves (300 mm long for unsize gloves)
- Our Packing room flooring is epoxy coated and all our packing tables are stainless steel ensuring supreme quality and a clean room atmosphere.



- We are an **ISO 9001** certified company. We have **EN46002** for medical device and comply with **EN455 Parts 1, 2 and 3**. All our products are registered with **FDA** including protein labeling claim. We adhere to the international standards and have state-of-the-art labs and equipment.



Products Information

- **Products currently offered: (Textured or Smooth Surface)**

- Powdered Examination Gloves
- Powder Free Chlorinated Examination Gloves
- Powder Free Polymer Coated Examination Gloves

- **Sizes**

- All sizes are available
 - XS, S, M, L, XL

- **Pricing Information**

Prices are based on FOB from Songkhla, Thailand per 1000 gloves.
(1 Carton = 10 boxes packed with 100 gloves)

Powdered Examination Gloves AQL 1.5 @ - **Ask for Quote**
Powder Free Chlorinated Examination Gloves - **Ask for Quote**

Powder Free Chlorinated Examination Gloves
with Low Protein Labeling - **Ask for Quote**

Polymer Coated Powder Free Examination Gloves - **Ask for Quote**



Contact Info:

Marketing Department
500 Moo4, T. Kuanlang,
Hatyai, Songkhla 90110
Tel: (66 74)474730-2
Fax: (66 74)474734
Mobile: (661) 839-3922
E-mail: inquiry@InnovativeGloves.net
Website: www.InnovativeGloves.net

We are big enough to Provide and small enough to Care !