

AVI NITRI SKIN PRO PRODUCT SPECIFICATIONS

Name of Product	3 Mil Powder free Nitrile Exam Gloves.
Name of Manufacturer	Innovative Gloves Co. Ltd. Hatyai, Thailand.
Base Material	Nitrile.
Characteristics	Non-Sterile, Ambidextrous Beaded Cuff.
Available Sizes	X-Small, Small, Medium, Large, X-Large, XX-Large.
Surface	Finger Textured.
Surface Treatment	Polymer Coating.
Color	White/Blue/Black/Violet.
Total Residue Content	< 2.0 Milligrams per glove.



Physical Specifications:

Size	Length(mm/inch)	Thickness(mm) Single Wall			Width(mm)
		Tip	Palm	Cuff	
X-Small	240±5	0.09±0.01	0.06±0.01	0.05±0.01	75±5
Small	240±5	0.09±0.01	0.06±0.01	0.05±0.01	85±5
Medium	240±5	0.09±0.01	0.06±0.01	0.05±0.01	95±5
Large	240±5	0.09±0.01	0.06±0.01	0.05±0.01	105±5
X-Large	240±5	0.09±0.01	0.06±0.01	0.05±0.01	115±5
XX-Large	240±5	0.09±0.01	0.06±0.01	0.05±0.01	120±3

Tensile (ASTM D412 & EN 455-2)

Before Aging

Tensile(MPA)	Elongation(%)	Force At Break(N)
22 Mpa Min.	500 Min.	≥ 6 N

After Accelerated Aging

Tensile(MPA)	Elongation(%)	Force At Break(N)
22 Mpa Min.	400 Min.	≥ 6 N




Packaging

Packaging (Normal)	100±2 pieces in one Dispenser Box,10 Dispenser Box in one Carton.
Labeling on Dispenser (Normal)	Product Name, Size, Lot no. Quantity, Manufacturer's address.
Labeling on Carton (Normal)	Product Name, Size, Lot no. Quantity, Country of Origin.

Container Loadability (No MOQ in manufacturer owned brand in standard Packaging)

20 FCL=2150 for 1000x10 Packing. 40 FCL=1360 for 1000x10 Packing.

Certification



EN420:2003+A1:2009



EN374-2:2003