


BLACK BEAUTY PRODUCT SPECIFICATIONS

Name of Product	Powder free Latex Gloves.
Name of Manufacturer	Innovative Gloves Co. Ltd. Hatyai, Thailand.
Base Material	Natural Rubber Latex.
Characteristics	Non-Sterile, Ambidextrous Beaded Cuff.
Available Sizes	X-Small, Small, Medium, Large, X-Large, XX-Large.
Surface	Smooth.
Surface Treatment	Mild Washing with Chlorine Solution and Extensive Washing.
Color	Black.
Total Residue Content	< 2.0 Milligrams per glove.
Total Extractable Protein content	< 50 μ gm/gm of glove.



Physical Specifications:

Size	Length(mm/inch)	Thickness(mm) Single Wall			Width(mm)
		Tip	Palm	Cuff	
X-Small	240 \pm 5	0.14 \pm 0.01	0.11 \pm 0.01	0.09 \pm 0.01	75 \pm 5
Small	240 \pm 5	0.14 \pm 0.01	0.11 \pm 0.01	0.09 \pm 0.01	85 \pm 5
Medium	240 \pm 5	0.14 \pm 0.01	0.11 \pm 0.01	0.09 \pm 0.01	95 \pm 5
Large	240 \pm 5	0.14 \pm 0.01	0.11 \pm 0.01	0.09 \pm 0.01	105 \pm 5
X-Large	240 \pm 5	0.14 \pm 0.01	0.11 \pm 0.01	0.09 \pm 0.01	115 \pm 5
XX-Large	240 \pm 5	0.14 \pm 0.01	0.11 \pm 0.01	0.09 \pm 0.01	120 \pm 3

Tensile (ASTM D412 & EN 455-2)

Before Aging

Tensile(MPA)	Elongation(%)	Force At Break(N)
21 Mpa Min.	700 Min.	≥ 6 N

After Accelerated Aging

Tensile(MPA)	Elongation(%)	Force At Break(N)
18 Mpa Min.	500 Min.	≥ 6 N

Packaging

Packaging (Normal)	100±2 pieces in one Dispenser Box,10 Dispenser Box in one Carton.
Labeling on Dispenser (Normal)	Product Name, Size, Lot no. Quantity, Manufacturer's address, Warning about Latex Allergy.
Labeling on Carton (Normal)	Product Name, Size, Lot no. Quantity, Country of Origin.

Container Loadability (No MOQ in manufacturer owned brand in standard Packaging)

20 FCL=1360 Cartons.	40 FCL=2740 Cartons.
----------------------	----------------------

Certification

