

ORANGE GRIPPER PRODUCT SPECIFICATIONS

Name of Product	7 Mil Powder free Nitrile Gloves Gripper.
Name of Manufacturer	Innovative Gloves Co. Ltd. Hatyai, Thailand.
Base Material	Nitrile.
Characteristics	Non-Sterile, Ambidextrous Beaded Cuff.
Available Sizes	Small, Medium, Large, X-Large, XX-Large.
Surface	Gripper Textured.
Surface Treatment	Polymer Coated.
Color	Orange.
Total Residue Content	< 2.0 Milligrams per glove.



Physical Specifications:

Size	Length(mm/inch)	Thickness(mm) Single Wall				Width(mm)
		Tip Smooth Surface	Tip Textured Surface	Palm Textured Surface	Cuff Smooth Surface	
Small	240±5	0.19±0.01	0.29±0.02	0.20±0.02	0.12±0.01	85±5
Medium	240±5	0.19±0.01	0.29±0.02	0.20±0.02	0.12±0.01	95±5
Large	240±5	0.19±0.01	0.29±0.02	0.20±0.02	0.12±0.01	105±5
X-Large	240±5	0.19±0.01	0.29±0.02	0.20±0.02	0.12±0.01	115±5
XX-Large	240±5	0.19±0.01	0.29±0.02	0.20±0.02	0.12±0.01	120±3

Tensile (ASTM D412 & EN 455-2)

Before Aging

Tensile(MPA)	Elongation(%)	Force At Break(N)
22 Mpa Min.	500 Min.	≥ 6 N

After Accelerated Aging

Tensile(MPA)	Elongation(%)	Force At Break(N)
22 Mpa Min.	400 Min.	≥ 6 N

Packaging

Packaging (Normal)	100±2 pieces in one Dispenser Box, 10 Dispenser Box in one Carton.
Labeling on Dispenser (Normal)	Product Name, Size, Lot no. Quantity, Manufacturer's address.
Labeling on Carton (Normal)	Product Name, Size, Lot no. Quantity, Country of Origin.

Container Loadability (No MOQ in manufacturer owned brand in standard Packaging)

20 FCL=1135 Cartons.	40 FCL=2740 Cartons.
----------------------	----------------------

Certification

